

Reference : CHANGE NOTICE No.: CN/15/06

Subject : Improvement in the identification of **mdi** products

Scope : This change notification will affect the following product categories supplied either non sterile or pre-sterilized:

- **mdi** Capsule Filters (All types and all sizes)
- **mdi** Cartridge Filters (All types and all sizes)
- **mdi** Centrifugal Filters (All types)
- **mdi** Funnels (All types)
- **mdi** Inline Filters (All types and all sizes)
- **mdi** Inline Vent Filters (All types and all sizes)
- **mdi** Membrane Disc Filters (All types and all sizes)
- **mdi** Stericheck: Closed Sterility Test System (All types)
- **mdi** Syringe Filters (All types and all sizes)
- **mdi** Vacufile Filters (All types)

Background:

In our efforts to provide the user better products with better identification and traceability, certain changes have been made in the identification of **mdi** products as mentioned above.

Need for Change:

To be able to inform the month of manufacture of products to our customers, labeling as well as printing of the products has been improved/modified to incorporate the month of manufacture.

How Does It Affect the User:

From a practical point of view the changes done will provide additional information to the user. All performance specifications and product attributes remain unchanged.

- The month of the manufacture of the product has been incorporated by suffixing an alphabet in the product lot number. Product identification/labeling as well as certificate of Quality will reflect the same. The example of proposed lot number are as below:

	Existing	Proposed
Product Category: Membrane Disc Filter (Pre-sterile)/Funnel		
Lot No.	1235	1235Z
Description	123 - Product batch number 5 - Last digit of year of manufacture	123 - Product batch number 5 - Last digit of year of manufacture Z - Month of manufacture*

	Existing	Proposed
Product Category: Capsule Filter/Cartridge Filter/Centrifugal Filter/Inline Filter/Inline Vent Filter/Stericheck/Syringe Filter/Vacufil Filter		
Lot No.	XY1245	XY1245Z
Description	X - Product Category Code Y - Membrane Code 123 - Product batch number 5 - Last digit of year of manufacture	X - Product Category Code Y - Membrane Code 123 - Product batch number 5 - Last digit of year of manufacture Z - Month of manufacture*
Product Category: Membrane Disc Filter (non sterile)		
Lot No.	XY123675	XY123675Z
Description	X - Membrane Code Y - Pore Size Code 123 - Last 3 digits of 1 st batch number 67 - Last 2 digits of last batch number 5 - Last digit of year of manufacture	X - Membrane Code Y - Pore Size Code 123 - Last 3 digits of 1 st batch number 67 - Last 2 digits of last batch number 5 - Last digit of year of manufacture Z - Month of manufacture*

* The month of manufacturing is as below:

A - January	B - February	C - March	D - April	E - May	F - June
G - July	H - August	I - September	J - October	K - November	L - December

Effective Date:

mdi product types as mentioned above manufactured after September 15, 2015 will have the said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



Quality Assurance
Advanced Microdevices Pvt. Ltd.