

Change Notice No.	Notification date	Implementation date
CN/19/010	July 01, 2019	August 01, 2019

Subject : Colour coding for mdi AseptiVac Vacuum Filters

Scope : This change notification will affect the products with following catalog numbers: (# can be any alphabet or numeral):

- AN#####
- AK#####

Background:

In our efforts to provide the user, products with better identification, certain improvements have been done in product identification of **mdi** AseptiVac Vacuum Filters with catalog numbers as mentioned above.

Need for Change:

Currently, all types of **mdi** AseptiVac Vacuum Filters with catalog numbers as mentioned above are manufactured using red colour caps. The different types of AseptiVac Vacuum Filters are identified with product related information such as product type, pore size, lot number etc. printed on product package and label pasted on product packaging box. It was observed that after opening the product package, it was difficult for the user to clearly identify different types of **mdi** AseptiVac Vacuum Filters. To make the product identification more user friendly, it has been decided to assign different colours to the caps used in different types of **mdi** AseptiVac Vacuum Filters, without changing the materials of construction. No change has been done in any dimensions or materials of construction or manufacturing process of **mdi** AseptiVac Vacuum Filters.

How Does It Affect the User:

From the user point of view, the changes done will improve the product identification under actual use. All performance specifications, product attributes and traceability will remain the same. Only the colour of the cap used in **mdi** AseptiVac Vacuum Filters is now re-defined as mentioned below:

- AseptiVac NS (Pore size: 0.2µm): Blue
- AseptiVac NS(Pore size: 0.45µm): Green
- AseptiVac KS (Pore size: 0.1µm, 0.2µm and 0.45µm): Red

Change Notice No.	Notification date	Implementation date
CN/19/010	July 01, 2019	August 01, 2019

Implementation of Change:

mdi AseptiVac Vacuum Filters with catalog numbers as mentioned above, manufactured after the implementation date will have the said changes. However, the available stock of **mdi** AseptiVac Vacuum Filters with only red colour cap will be received by you till the stocks last.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



Head - Quality Assurance
Advanced Microdevices Pvt. Ltd.