

Change Notice No.	Notification date	Implementation date
CN/A/23/005	September 22, 2023	October 22, 2023


Subject: Improvement of mdi SS Trolley and liner drum.

Scope: This change will affect mdi SS Trolley and liner drum with following catalog numbers: (X can be any alphabet or numeral):


- TSSXIXXXXXXXXXX01
- TSSXKXXXXXXXXXX01
- PEDXIXXXXXXXXXX01
- PEDXKXXXXXXXXXX01

Background:


To provide a better fitment for liner bags, liner drums and SS trolleys need to be more robust.



Old Design of SS Trolley



Proposed Design of SS Trolley



Detail:A:-
Bolts added on the bottom side to assemble wheel.

Change Notice No.	Notification date	Implementation date
CN/A/23/005	September 22, 2023	October 22, 2023

Need for Change:

The change will be done to improve the robustness of liner drum and SS trolleys.

The change in liner drum, will provide better fitment of liner bag in the drum and make it more convenient .

In SS trolleys, the detachable wheels will be used instead of welded wheels of trolley to make it user friendly.

How Does It Affect the User:

From a user point of view, the change will improve the robustness of liner drums and SS trolleys and in no way will affect the performance characteristic of the product under actual conditions of use. All performance specifications (except minor changes in dimensions), product attributes, material of construction, and traceability shall remain the same.

Effective Date:

The change in liner drum and SS trolley shall be implemented from October 22, 2023. However the available stock of the mdi SS trolley and liner drum will be shipped till stocks last.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



Head - Quality Assurance (Single Use System)
Advanced Microdevices Pvt. Ltd.