

| Change Notice No. | Notification date | Implementation date |
|-------------------|-------------------|---------------------|
| CN/A/23/006 | December 30, 2023 | January 30, 2024 |

Subject: Usage of Oetiker Clamps instead of wire ties in mdi Single Use Filter Assemblies and Single Use Bag Assemblies having tubing with Inner Diameter ≥ 12.7 mm.

Scope: This change notification will affect the Single Use Assemblies having tubing with Inner Diameter ≥ 12.7 mm.

Background:

Use of wire ties to secure large ID tubings connections may have possibility of leakage. In case of large tubing connection to provide better fitment and reliable connection, oetiker clamps shall be used instead of wire ties.

Need for Change:

Requirement of stronger and more reliable connection for large ID tubing with inner diameter ≥ 12.7 mm.

How Does It Affect the User:

From a user point of view, the change shall provide better fitment and leak proof connection with the use of oetiker clamp instead of wire ties. The change does not affect the quality of the product. All performance specifications and product attributes will remains the same.

Effective Date:

mdi Single Use Assemblies and Systems, released after the implementation date (January 30, 2024) will have the said changes. All the affected approved drawings shall be updated before manufacturing. In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



**Head - Quality Assurance (Single Use Systems)
Advanced Microdevices Pvt. Ltd.**