

Change Notice No.	Notification date	Implementation date
CN/A/22/005	November 09, 2022	December 09, 2022

Subject: A protective sleeve is provided on sampling port of mdj Single Use Assemblies, ASESS Sampling Systems and Single Use Components.

Scope: This change notification will affect the products with following catalog numbers: (X can be any alphabet or numeral):

- AXXXXXXXXXXXXXXXXX
- SSXXXXXXXXXXXXXX

Background:


The sampling port is an integral part of the mdj sampling systems as it provides a sterile pathway to draw samples by piercing the sampling port septum. A protective sleeve is provided on the septum of the sampling port.

Need for Change:

The protective sleeve will provide protection from any damage to the sampling port.

How Does It Affect the User:

From the user point of view, the change shall result in providing protection from any damage to the sampling port without impacting the quality of the product. All performance specifications, product attributes, material of construction, and traceability shall remain the same.



Change Notice No.	Notification date	Implementation date
CN/A/22/005	November 09, 2022	December 09, 2022

Effective Date:

mdi sampling system, manufactured after the implementation date (December 09, 2022) will have the said changes. However the available stock of the **mdi** Single Use Assemblies and Sampling Systems without protective sleeve will be shipped till stocks last.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



Head - Quality Assurance (Single Use Systems)
Advanced Microdevices Pvt. Ltd.