

Reference : CHANGE NOTICE No. CP/12/01

Subject : Modification in Disc sealing process of 1" Capsule Filter type DKL

Scope : This change notification will affect the following products:

- **mdi** capsule filter type DKL of 1" length (170 cm²)

Background:

In our efforts to provide the user better products with higher reliability, certain changes have been made in the Disc sealing process of **mdi** capsule filter type DKL of 1" length (170 cm²).

Need for Change:

To improve the productivity and reliability of production process of the capsule filters.

How Does It Affect the User:

From a user point of view the changes done only affect the visual appearance of the mini cartridge height inside the capsule body. It appears of lower height than the one with earlier process.

There is no change in mini cartridge design or membrane area.

- **Disc sealing Process** – Jointing of mini cartridge to the capsule cap is now done by direct sealing of mini cartridge in to the capsule cap instead of through an end disc earlier. This reduces the number of joints resulting in increased reliability and improved productivity.

Effective Date:

mdi Capsule filters type DKL of 1" length (170 cm²) manufactured after July 25, 2011 will be with said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'mdi@vsnl.com'.



Quality Assurance
Advanced Microdevices Pvt. Ltd.