

Reference : CHANGE NOTICE No. CP/13/09

Subject : Modification in identification on Capsule Filter of length 2", 5" and 8" (All types)

Scope : This change notification will affect the following products:

- All types of capsule filters of length 2", 5" and 8"

Background:

In our efforts to provide the user better products with better identification and traceability, certain changes have been made in the identification on **mdi** capsule filter (all types) of length 2", 5" and 8".

Need for Change:

To improve the identification and traceability of the capsule filters under different conditions of use.

How Does It Affect the User:

From a user point of view the changes done will improve the identification, traceability and ease-of-use under actual conditions of use. All performance specifications, product attributes and main content of the identification remain the same.

- **Identification on capsule filter** – Laser etching method is now used to print the identification details on capsule filter which was earlier done by embossing on the capsule filters. This will improve readability under different conditions of use.

OLD (Embossing)



NEW (Laser Etching)



Effective Date:

mdi Capsule Filters (all types) of length 2", 5" and 8" manufactured after November 15, 2013 will be with said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'mdi@vsnl.com'.



Quality Assurance

Advanced Microdevices Pvt. Ltd.