

**Reference :** CHANGE NOTICE No. MB/11/04

**Subject :** Change in mdi Closed Sterility Test Device type Stericheck-SVP1 and Stericheck-SVP3 (All membrane types)

**Scope :** This change notification will affect the following component of mdi Stericheck:

- Needle - New improved design for liquid needle

**Background:**

In our efforts to provide the user better products with more convenience, certain changes have been made in the needle of mdi Closed Sterility Test Device type **Stericheck-SVP1 and Stericheck-SVP3** (All membrane types).


**Need for Change:**

As above.

**How Does It Affect the User:**

From a practical point of view this change will improve the needle performance under actual use. All performance specifications and product attributes remain unchanged.

- **Stericheck needle for SVP1 and SVP3:** A new improved design of liquid needle has been made which will eliminate the cutting of rubber bung during actual use.



**Effective Date:**

mdi Closed Sterility Test Device type **Stericheck-SVP1 and Stericheck-SVP3** (all membrane types) manufactured after April 15, 2011 will have said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'mdi@vsnl.com'.

**Quality Assurance**

**Advanced Microdevices Pvt. Ltd.**