

Reference : CHANGE NOTICE No. : MB/12/01

Subject : Rationalization of pore size of Membrane Disc Filters type CN/0.22 µm (all diameters)

Scope : This change notification will affect the following product:

- **mdi** 0.22 µm Membrane Disc Filter Type CN (All diameters)

Background:

In our efforts to rationalize our products, the pore size of 0.22 µm Membrane Disc Filters Type CN (all diameters) has been designated as 0.2 µm.

Need for Change:

As above.

How Does It Affect the User:

From a practical point of view this change will only change the pore size designation of membrane disc filters. All performance specifications and product attributes including bacterial retention abilities remain unchanged.

- The rationalized pore size of membrane disc filter type CN (all diameters) is now 0.2 µm instead of 0.22 µm earlier.

Effective Date:

mdi 0.22 µm membrane disc filters type CN (all diameters) shipped/ordered/processed after May 28, 2012 will be with said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'mdi@vsnl.com'.



Quality Assurance
Advanced Microdevices Pvt. Ltd.