

Reference : CHANGE NOTICE No.: CN/15/05

Subject : Improved design of end connection type 1/4" - 3/8" Hose barb (Code - B) used with **mdi** Inline Filters and Inline Vent Filter

Scope : This change notification will affect the following products:  
**mdi** 50mm Inline Filters and Inline Vent Filters with end connection type 1/4" - 3/8" Hose barb (Code - B) - all pore sizes

**Background:**


In our efforts to provide the user better products, certain changes have been made in the design of end connection type 1/4" - 3/8" Hose barb (Code - B) used with **mdi** Inline Filters and Inline Vent Filter as mentioned above.

**Need for Change:**

As above.

**How Does It Affect the User:**

From a practical point of view the changes done will improve the fitment of respective tubings under actual use including fitment effectiveness during autoclaving. All performance specifications and product attributes remain unchanged, except the improved fitment.



**Effective Date:**

**mdi** Inline Filters and Inline Vent Filters (as mentioned above) manufactured after August 08, 2015 will have said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



Quality Assurance  
Advanced Microdevices Pvt. Ltd.