

Change Notice No.	Notification date	Implementation date
CN/16/11	05 ^h August 2016	05 th September 2016

Subject : Change in Lot numbering system of **mdi** pre-sterilized Syringe Filters (All types)

Scope : This change notification will affect the lot number of products with following catalog numbers: (* can be any alphabet or numeral):

- SMKG*****2**
- SYCA*****2**
- SYKG*****2**
- SYN*****2**
- SYNN*****2**
- SYNU*****2**
- SYPG*****2**
- SYPL*****2**
- SYPL*****3**
- SYPN*****2**
- SYPP*****2**
- SYTF*****2**
- SYTG*****2**
- SYVF*****2**
- SYVF*****3**
- SYVG*****2**

Background:

In keeping with our continuous efforts to provide the user products with easy and clear identification, certain changes have been made in the lot number system of **mdi** pre-sterilized Syringe Filters with catalog numbers as mentioned above.

Need for Change:

Lot number identification on packaging label of **mdi** pre-sterilized Syringe Filters includes product code, membrane code, last 3 digits of production batch number, last digit of manufacturing year, manufacturing month code, pouch sealing sheet number and sterilization cycle number. Pouch sealing sheet number and sterilization cycle number are three digits each. The possibility of interchange in these two values could not be ruled out resulting in wrong traceability. The pouch sealing sheet number is a production step which can be traced back through production batch number. Hence, it is not necessary to include it in the Lot Number. Thus the lot number identification on packaging label of **mdi** pre-sterilized Syringe Filters has been changed.



Change Notice No.	Notification date	Implementation date
CN/16/11	05 th August 2016	05 th September 2016

How Does It Affect the User:

From a practical point of view, this change only affects the lot numbering system printed on packaging label. All performance specifications, product attributes, materials of construction, product packaging and traceability remain unchanged.

Lot number identification on packaging label has been changed to lot number and sterilization cycle number (Ster. No.) only. Pouch sealing sheet number has been removed. Product packaging label will change accordingly (table-1). **mdi** pre-sterilized Syringe Filters manufactured after the implementation date as mentioned above will have the said changes.

Table - 1

	EXISTING	PROPOSED
Lot No.	Lot No.: SN5736E-123-001	Lot No.: SN5736E Ster. No. 001
Lot number description	'S' is the product code 'N' is the membrane code '573' are last 3 digits of production batch number '6' is the last digit of respective year 'E' is the code for month of manufacture '123' is the pouch sealing sheet No. '001' is sterilization cycle No.	'S' is the product code 'N' is the membrane code '573' are last 3 digits of production batch number '6' is the last digit of respective year 'E' is the code for month of manufacture '001' is sterilization cycle No.
Product packaging Label	 <p>mdi SYRINGE FILTER SINGLE USE ONLY Cat. No. : SYNNO601MNXX204 Type : SY25NN-S (Nylon) Pore Size : 0.2 µm Diameter : 25 mm Quantity : 100 Nos. Lot No. : SN5736E-123-001</p> <p>STERILE EO 2019 - 07</p> <p>ADVANCED MICRODEVICES PVT. LTD. 21, IND. AREA, AMBALA CANTT. - 133 006 INDIA www.mdimembrane.com</p>	 <p>mdi SYRINGE FILTER SINGLE USE ONLY Cat. No. : SYNNO601MNXX204 Type : SY25NN-S (Nylon) Pore Size : 0.2 µm Diameter : 25 mm Quantity : 100 Nos. Lot No. : SN5736E Ster. No. : 001</p> <p>STERILE EO 2019 - 07</p> <p>ADVANCED MICRODEVICES PVT. LTD. 21, IND. AREA, AMBALA CANTT. - 133 006 INDIA www.mdimembrane.com</p>

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.

