

Change Notice No.	Notification date	Implementation date
CN/17/02	10 th January 2017	10 th February 2017

Subject : Improved design of of mdi Funnel (Funnel Base)

Scope : This change notification will affect the products with following catalog numbers (where * can be any alphabet or numeral):

- FMCN*****
- FMCP*****
- FMEN*****
- FMMC*****
- FMNN*****
- FMVF*****

Background:

In our efforts to improve the ease of use of the product, certain changes have been made in the design of Funnel base of **mdi** Funnel with catalog numbers as mentioned above.


Need for Change:

During use of **mdi** Funnel, the membrane is to be cultured with growth medium for incubation. To perform this process step, the membrane assembled in funnel base has to be lifted and transferred to respective growth medium plates to facilitate the growth of microorganism retained on the membrane. It was observed that It is sometimes difficult to lift the membrane from funnel base. To address this observation, design of funnel base of **mdi** Funnel has been improved by providing a groove for easy lifting of membrane. No change has been done in any critical dimensions of the Funnel. The new design of **mdi** Funnel has been verified for critical product performance parameters including fitment with respective components.

How Does It Affect the User:

From a practical point of view the changes done will improve the ease of use of **mdi** Funnel. All performance specifications, product attributes, materials of construction and fitment with respective components will remain unchanged.

Change Notice No.	Notification date	Implementation date
CN/17/02	10 th January 2017	10 th February 2017



Implementation of Change:

mdi Funnel with catalog numbers as mentioned above, manufactured after the implementation date (**10th February 2017**) will have the said changes. However, the available stock of **mdi** Funnel with old design of Funnel base will be received by you till the stocks last.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



Head - Quality Assurance
Advanced Microdevices Pvt. Ltd.