

Change Notice No.	Notification date	Implementation date
CN/17/21	November 01, 2017	December 01, 2017

Subject : Addition of Microbial Recovery Test in lot release criteria in Certificate of Quality (CoQ) of mdi Funnel

Scope : This change notification will affect the products with following catalog numbers: (# can be any alphabet or numeral):

- FMCN#####
- FMEN#####
- FMMC#####
- FMNN#####
- FMVF#####

Background:

In our efforts to provide the user latest and current information regarding the performance parameters of the products they are using, it has been decided to make certain addition in the contents of the Certificate of Quality (CoQ) of **mdi** Funnels with catalog numbers as mentioned above.

Need for Change:

mdi Funnels with catalog numbers as mentioned above are validated for Microbial Recovery as per ISO 7704-1985 (E): Evaluation of membrane filters used for microbiological analyses. To further enhance the verification of Microbial Recovery, it has been decided to perform Microbial Recovery of *Escherichia coli* (ATCC 8739) on every manufactured lot of **mdi** Funnels. To incorporate this updated information, the CoQ of **mdi** Funnels with catalog numbers as mentioned above has been updated with inclusion of Microbial Recovery of *Escherichia coli* (ATCC 8739) in lot release criteria.

Change Notice No.	Notification date	Implementation date
CN/17/21	November 01, 2017	December 01, 2017

How Does It Affect the User:

From a practical point of view this change only affects the CoQ of above mentioned product types by providing additional information in the form of Microbial Recovery of *Escherichia coli* (ATCC 8739) in lot release criteria. All product attributes and other performance specifications remain unchanged. The accompanying CoQ will reflect the said change.

Implementation of Change:

mdi Funnel with catalog numbers as mentioned above, manufactured after the implementation date (**December 01, 2017**) will have the said changes. However, the available stock of **mdi** Funnel with catalog numbers as mentioned above having CoQ without Microbial Recovery of *Escherichia coli* (ATCC 8739) will be received by you till the stocks last.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



Head - Quality Assurance
Advanced Microdevices Pvt. Ltd.