

Change Notice No.	Notification date	Implementation date
CN/19/003	March 01, 2019	April 01, 2019

Subject : Change in Funnel capacity of mdi Vacuum Filter Unit for Glass Bottles

Scope : This change notification will affect the products with following catalog numbers: (# can be any alphabet or numeral):

- AKX7####F##G###
- AKX7####Z##G###

Background:

In our efforts to provide convenience to the user/customer, certain changes have been made in the funnel of **mdi** Vacuum Filter Unit for glass bottles having catalog numbers as mentioned above. The capacity of the funnel used in above **mdi** Vacuum Filter Unit for glass bottles has been upgraded from existing 500 mL to 1000 mL.

Need for Change:

Currently, a 500 mL funnel is being supplied as component of **mdi** Vacuum Filter Unit for glass bottles having catalog numbers as mentioned above. The funnel is used to hold the unfiltered fluid and the required volume of fluid to be filtered is sometimes more than 500 mL. To address this issue, capacity of the funnel used in above **mdi** Vacuum Filter Unit for glass bottles has been increased from 500 mL to 1000 mL. No change has been done in materials of construction, manufacturing process, product attributes and product specifications.

How Does It Affect the User:

From a practical point of view this change only affect the funnel supplied with **mdi** Vacuum Filter Unit for glass bottles with catalog numbers as mentioned above by an increase in capacity to 1000 mL. All product attributes and performance specifications remain unchanged. However, product packaging will be changed accordingly.

Change Notice No.	Notification date	Implementation date
CN/19/003	March 01, 2019	April 01, 2019

Implementation of Change:

mdi Vacuum Filter Unit for glass bottles having catalog numbers as mentioned above, manufactured after the implementation date (**April 01, 2019**) will have the said changes. However, the available stock of **mdi** Vacuum Filter Unit for glass bottles with catalog numbers as mentioned above having funnel with capacity 500 mL will be received by you till the stocks last.

In case you have any queries, please feel free to contact our Technical Support Team at 'info@mdimembrane.com'.



Head - Quality Assurance
Advanced Microdevices Pvt. Ltd.