

Reference : CHANGE NOTICE No. CP/13/05

Subject : Improved design of Protective Caps for vent/Drain valve and Inlet/Outlet connection (½" Stepped Hose Barb only) used with mdi Capsule Filters

Scope : This change notification will affect the following products:

- **mdi** Gamma Sterilizable and Gamma Irradiated Capsule Filters (All types with sizes 2", 5" and 8")
 - Protective caps for vent/drain valve – All end connections
 - Protective caps for Inlet/Outlet connection – ½" Stepped Hose Barb only

Background:

In our efforts to provide the user better products, certain changes have been made in the design of Protective Cap for vent/Drain valve and Inlet/Outlet connection (½" Stepped Hose Barb only) used with **mdi** Capsule Filters types mentioned above.



Need for Change:

As above.

How Does It Affect the User:

From a practical point of view the changes done will improve the robustness of product during transportation. All performance specifications and product attributes remain unchanged.

- **Protective caps for vent/drain valve** - protective caps are now inserted in the vent/drain valve, which were earlier assembled over the vent/drain valve. This will reduce loosening of the protective caps during transportation.
- **Protective caps for ½" Stepped Hose Barb Inlet/Outlet connection** - protective caps are now inserted in the inlet/outlet connection, which were earlier assembled over the inlet/outlet connection. This will reduce loosening of the protective caps during transportation.

<p align="center"><u>OLD PROTECTIVE CAPS</u> (1/2" Stepped Hose Barb Inlet/outlet and vent/drain valve)</p>	<p align="center"><u>NEW PROTECTIVE CAPS</u> (1/2" Stepped Hose Barb Inlet/outlet and vent/drain valve)</p>
 <p>The image shows a white cylindrical membrane filter capsule with a red cap on top. It has two side ports with red caps and a bottom port with a red cap. A label on the side reads: 'mdi CAPSULE FILTER DKS-S 0.2µm Cat. No. DKS-S0200X000 Advanced Microdevices P. Ltd www.mdimembrane.com'.</p>	 <p>The image shows a white cylindrical membrane filter capsule with a red cap on top. It has two side ports with red caps and a bottom port with a red cap. A label on the side reads: 'mdi CAPSULE FILTER DKS-S 0.2µm Cat. No. DKS-S0200X000 Advanced Microdevices P. Ltd www.mdimembrane.com'.</p>

Effective Date:

mdi Capsule filters (types and sizes as mentioned above) manufactured/shipped after October 25, 2013 will have said changes.

In case you have any queries, please feel free to contact our Technical Support Team at 'mdi@vsnl.com'.



Quality Assurance
Advanced Microdevices Pvt. Ltd.