

# Certificate of Quality

The Polycarbonate Quick Connector has been manufactured in a **mdi** facility in compliance with **ISO 9001** regulations.

## Polycarbonate Quick Connector (Coupling Body)

Catalog No. : QCCX95CBXXXX104  
Tubing Size : 9.5 mm (3/8")  
Lot Number : PC0278K

### SPECIFICATION

Material of Construction : Polycarbonate  
Operating Pressure : Vacuum at 60 psi  
Operating Temperature : - 60 to 150 °C  
Sterilization : 5 autoclaving cycles at 121 °C for 30 minutes each

### LOT RELEASE CRITERIA

Each lot has been sampled, tested and released by the Quality Assurance Department for the following characteristics:

- o Visual inspection
- o Pressure test

### VALIDATED FOR

**Biosafety** : Passes Biological Reactivity Tests, *In Vivo* for Class VI plastic as described in USP <88>.  
**Cytotoxicity** : Passes Biological Reactivity Tests, *In Vitro* as described in USP <87>.  
**Fiber Release** : Complies with USFDA 21CFR 211.72.

### CUSTOMER SUPPORT

**mdi** offers its unique interdisciplinary skills to provide solutions to specific problems. Please contact our factory or the local application specialist.



Head of Quality Assurance

Issue Date: 29-Nov-18

**Advanced Microdevices Pvt. Ltd.**  
21, Industrial Area, Ambala Cantt, INDIA,  
Tel: +91-171-2699290/ 2699274  
Website: [www.mdimembrane.com](http://www.mdimembrane.com)  
Email: [info@mdimembrane.com](mailto:info@mdimembrane.com)

**An ISO 9001 Company**